



German Standard Medical Gas Outlets



We are a company somewhere in China, we have our own factory, all products are factory direct sales.

DFTT® German Standard Medical Gas Outlets is one of our most popular products, which have been sold all over the world, such as Eastern Europe, Southeast Asia, America and so on. German Standard Gas Outlets are of high quality and low price, each of them uses superior materials and advanced technology, so they are very popular with customers. We hope that we can also build a good relationship between us, whether long-term or short-term, as long as you can trust us, we will make you satisfied.

Product Description

DFTT® German Standard Medical Gas Outlets consist of copper material, waterproof and fireproof in daily life, not endangering life safety, so that patients rest assured. German Standard Gas Outlets are irregular in shape, but this does not interfere with normal operation. Generally speaking, German Standard Gas Outlets are equipped with bed Head unit to control oxygen flow and ensure smooth breathing.

Product Parameter

Product name	Shape	material	size	color
German Standard Gas Outlets	irregular shape	Cooper	7*6*6CM	Golden

German Standard Gas Outlets Feature And Application



German Standard Medical Gas Outlets have a variety of shapes, due to the diversity of Gas inlet leads to a variety of connections, the use of color to distinguish the Gas inlet, easy to connect the Gas pipeline, special gas, greatly ensure the safety, avoid tragedy. German Standard Gas Outlets are commonly used for connecting the bed head unit to the inhaler in hospital wards. The German Standard Gas Outlets have built-in access valves for convenient and simple maintenance without disassembly of the panel. In addition, the direction of import and export can be adjusted to achieve real time and labor saving.

German Standard Gas Outlets Details







German standard gas outlets DMG09-3

The biggest advantage of DFTT® German Standard Medical Gas Outlets is that they can be pressed with one hand and do not require great effort, which is impossible with other similar products. From the details, the surface after fine decoration treatment, simple and beautiful, small and portable. German Standard Gas Outlets are suitable for German plugs, preventing other plugs from being inserted by mistake, thus reducing safety hazards. German Standard Gas Outlets through tens of thousands of plugging experiments and air tightness tests, safe and reliable, meet the national quality standards and passed the CE certification.